

## AMPUTATED AND AVULSED PARTS

1. Adult/Pediatric Routine Trauma Care
2. If gross contamination, gently rinse stump and amputated part with NORMAL SALINE.
3. Control bleeding with direct pressure. Utilize hemostatic gauze if available. If needed, apply a tourniquet.
4. Cover stump with sterile NS moistened dressing, additional sterile dressing, and wrap with elastic bandage (ACE wrap) providing uniform pressure over entire stump.
5. Amputated part:
  - A. In case of incomplete amputation, splint in physiologic position.
  - B. Wrap part in NS moistened gauze and place in a plastic bag/container.
  - C. Place bag/container in in a larger bag or container with cold packs or ice and water.
  - D. Do not allow parts to come in direct contact with ice/cold packs or submerge in NS.